



PostureRay® Minimum Requirements

The following computer system and software specifications are the minimal requirements for the PostureRay® X-Ray software. Most new computers purchased within the last two years meet these minimal requirements, which include:

Minimum Stand-Alone PostureRay® Hardware Requirements:

Component	Specification
CPU	Intel Core i3 to i7 Recommended
RAM	>4 GB Minimum Required
Storage	>100 GB Available for Images and Reports
Video	>17"-24" Monitor at 1024 x 768 Resolution
Optical Drive	CD/DVD for Referring Physician Images
Other	Internet Connection – Broadband Cable Recommended

Minimum PostureRay® Software Requirements:

Component	Specification
Operating System*	Minimum Windows 7 Service Pack 1 Required. Professional Version for Windows 8 and 10 with Current Patching Recommended.
PDF Reader	Adobe Reader DC https://get.adobe.com/reader/
Internet Browser	Internet Explorer 11 or later, Chrome, or FireFox
Firewall	Network or Windows Firewalls must allow PostureRay to contact the Internet
Other	Windows .NET Framework 4.5

* MAC users must have Parallels Desktop, Bootcamp or VMWare Fusion installed with a minimum of Windows 7 SP1 with all available updates installed.

Minimum PostureRay® Network Requirements (not sold nor maintained through PostureCo, Inc):

Component	Specification
Server	See Hardware/Software Specifications Above. For a Stand-Alone Server, Windows 2008 Minimum is Required, Windows Server 2012 or 2016 Recommended
Marking Station/Viewing Station	See Specifications Above. Can be PC, Laptop, or Tablet Computer
Network	Standard TCP/IP Networking. NOTE – Wired Networks Preferred for Stability and Speed
Remote Connectivity	NOTE - Terminal Services is not Supported!

Other optional Equipment not sold or maintained through PostureCo, Inc.

Component	Specification
Printer	Color Printer is Recommended
DICOM x-ray Scanner*	Scan X-Rays into .DCM Format. DICOMS should be 12-16 bit High Quality JPEG2000 Compressed
Digital Camera	3 Megapixel or Better Digital Camera – NOT PHONE CAMERA
X-Ray Digitized Machine*	Connectivity to a CR and/or DR Imaging System that uses PACS. Important Message: PostureRay is not a PACS Receiving Software

* Images must be *.JPG, *.BMP, or true DICOM with *.DCM extensions

PostureCo, Inc.

3152 Little Road, Suite 161, Trinity, FL 34655

Tele: 866-577-7297 Fax 727-369-8810 Email: info@postureco.com